NovaLisa™

**Liquor diagnostics
Borrelia burgdorferi**

IgG/ IgM

CE

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Page
Seite

2 to 6
7 bis 12

Product Number: BORL0040

1. GENERAL ASPECTS

In case of clinically suspected neuroborreliosis, the detection of intrathecally produced *Borrelia*-specific IgG and IgM antibodies is indicated. For this purpose serum-CSF pairs are analysed in parallel in the quantifiable NovaLisa™ *Borrelia burgdorferi* IgG/ IgM ELISA. The result can then be confirmed by Blot.

First, serum and CSF are included in the standard *Borrelia* screening. If the CSF is positive, the values for total IgG and/ or IgM and albumin must be requested for serum and CSF as the basis for the detection of the *Borrelia* specific antibody index.

By this method based on the *Borrelia*-specific antibody index, an intrathecally produced antibody fraction can be differentiated from a fraction that has passed the blood-CSF barrier passively to enter the subarachnoid space. Applying the limiting value QLim from the quotient diagram according to Reiber ensures that no false-negative antibody index values will be calculated in case of a polyspecific intrathecally immune reaction.

In cases of serum-CSF pairs showing a pathologically elevated antibody index in the screening test, confirmation with Blot is recommended.

2. INTENDED USE

The NovaLisa™ *Borrelia burgdorferi* IgG and IgM ELISA are intended for the quantitative detection of intrathecally produced human IgG and IgM antibodies against *Borrelia burgdorferi*.

3. SPECIMEN COLLECTION AND PREPARATION

Serum and CSF must be withdrawn from the patient at the same day. Parallel testing of serum and CSF in the same run is recommended. To increase diagnostic safety, patient sera and CSF should be tested in duplicate.

3.1. Sample Dilution

IgG

Serum 1:101 The controls of the ELISA kit are ready to use. Patient samples have to be diluted 1:101 with the yellow sample diluent. Mix well! 100 µl of each (ready to use controls and diluted samples) are pipetted per well.

CSF 1:4 75 µl of the yellow sample diluent are pipetted directly into the well. Then 25 µl CSF are added. Mix well! Please be aware of the correct order.

IgM

Serum 1:101 The controls of the ELISA kit are ready to use. Patient samples have to be diluted 1:101 with the green sample diluent. Mix well! 100 µl of each (ready to use controls and diluted samples) are pipetted per well.

CSF 1:2 50µl of the green sample diluent are pipetted directly into the well. Then 50 µl CSF are added. Mix well! Please be aware of the correct order.

4. ASSAY PROCEDURE

The further test procedure is performed according to the instruction for use, NovaLisa™ *Borrelia burgdorferi* IgG and/ or IgM.

5. EVALUATION

The test is valid if the criteria for the control sera are fulfilled according to the instructions for use of the ELISA NovaLisa™ *Borrelia burgdorferi* IgG and/ or IgM.

Calculation of threshold values for sera and CSF for the qualitative evaluation:

IgG/ IgM serum: threshold value = OD cut-off control

IgG/ IgM CSF: threshold value = 0.5 x OD cut-off control

5.1. Interpretation of results

IgG result	Interpretation and further actions
OD _{Serum} 1:101 < threshold value serum and OD _{CSF} 1:4 < threshold value CSF	No specific IgG antibodies in serum detected. No specific IgG antibodies in CSF detected. No indication of neuroborreliosis No determination of antibody index (AI) IgG
OD _{Serum} 1:101 ≥ threshold value serum and OD _{CSF} 1:4 < threshold value CSF	Detection of specific IgG antibodies in serum No specific IgG antibodies in CSF detected. Indication of an advanced or past Borrelia infection No indication of neuroborreliosis No determination of antibody index (AI) IgG
OD _{Serum} 1:101 < threshold value serum und OD _{CSF} 1:4 ≥ threshold value CSF	No specific IgG antibodies in serum detected. Detection of specific IgG antibodies in CSF Neuroborreliosis is possible Determination of antibody index (AI) IgG
OD _{Serum} 1:101 ≥ threshold value serum und OD _{CSF} 1:4 ≥ threshold value CSF	Detection of specific IgG antibodies in serum Detection of specific IgG antibodies in CSF Neuroborreliosis is possible Determination of antibody index (AI) IgG

IgM results	Interpretation and further actions
OD _{Serum} 1:101 < threshold value serum and OD _{CSF} 1:2 < threshold value CSF	No specific IgM antibodies in serum detected. No specific IgM antibodies in CSF detected. No indication of neuroborreliosis No determination of antibody index (AI) IgM
OD _{Serum} 1:101 ≥ threshold value serum and OD _{CSF} 1:2 < threshold value CSF	Detection of specific IgM antibodies in serum No specific IgM antibodies in CSF detected. Indication of a fresh Borrelia infection No indication of neuroborreliosis No determination of antibody index (AI) IgM
OD _{Serum} 1:101 < threshold value serum und OD _{CSF} 1:2 ≥ threshold value CSF	No specific IgM antibodies in serum detected. Detection of specific IgM antibodies in CSF Neuroborreliosis is possible Determination of antibody index (AI) IgM
OD _{Serum} 1:101 ≥ threshold value serum und OD _{CSF} 1:2 ≥ threshold value CSF	Detection of specific IgM antibodies in serum Detection of specific IgM antibodies in CSF Neuroborreliosis is possible Determination of antibody index (AI) IgM

The test can be terminated if serum and CSF give negative IgG and IgM results in screening. The Borrelia specific antibody index cannot be determined.

In case the CSF OD exceeds the threshold value, the procedure continues. In the following , IgG and IgM tests will be described in parallel. However, the antibody index is to be determined only for the immunoglobulin class in which a positive screening result for the CSF has been obtained.

5.2. Calculation

5.2.1. Calculation of the specific IgG resp. IgM in NTU/ml

For the conversion of the ODs of CSF and serum into units, extinction values ≤ 2.500 are required.

If the OD value for the serum or CSF is larger than 2.500 or even outside the measurable range, the serum-CSF pair must be retested together in another run after appropriately higher dilution. The following dilutions are recommended:

IgG

OD	Serum dilution	CSF dilution	
		Preparation in the well	
Standard dilution	1:101	1:4	75 µl SD* + 25 µl CSF

Result of standard dilution (OD):	Recommendation for further dilutions:		
OD ≤ 2,500	No further dilution		
2,500 < OD ≤ 3,000	1:500	1:10 1:20	90µl SD* + 10 µl CSF 95 µl SD* + 5 µl CSF
OD > 3,000	1:1000; 1:2000; 1:4000	Dilution in a polypropylen (PP) tube 1:40; 1:200; 1:1000	

IgM

OD	Serum dilution	CSF dilution	
		Preparation in the well	
Standard dilution	1:101	1:2	50 µl SD* + 50 µl CSF

Result of standard dilution (OD):	Recommendation for further dilutions:		
OD ≤ 2,500	No further dilution		
2,500 < OD ≤ 3,000	1:500	1:5	80µl SD* + 20 µl CSF
OD > 3,000	1:1000; 1:2000; 1:4000	Dilution in a polypropylen (PP) tube 1:10; 1:20; 1:40	

* SD = sample diluent

The specific IgG resp IgM has to be calculated from OD values of serum and CSF ≤ 2.500.
If several values are in this range, the value with the OD closest to 1.000 is selected.

Note: As the cut-off in this formula the value of the cut-off control (cut-off for serum) of the corresponding IgG or IgM determination is to be used for serum and for CSF. The threshold value for CSF from the qualitative evaluation is not valid here.

$$(1) \quad \text{spec. IgG resp. IgM serum resp. CSF [NTU/ml]} = \frac{\text{OD}_{\text{Serum resp. Liquor}} \times 10 \times \text{dilution factor}}{\text{cut-off} \times 101}$$

5.2.2. Calculation of the specific IgG- resp. IgM - quotient (Qspec IgG/ Qspec IgM)

- (2) $Q_{\text{spec}}(\text{IgG}) = \text{spec. IgG CSF [NTU/ml]} / \text{spec. IgG serum [NTU/ml]}$
- (3) $Q_{\text{spec}}(\text{IgM}) = \text{spec. IgM CSF [NTU/ml]} / \text{spec. IgM serum [NTU/ml]}$

5.2.3. Calculation of the quotients QIgG, QIgM, QAlb

Total IgG, total IgM and albumin in both serum and CSF are determined by using a suitable method.
From the values obtained, the IgG quotient Q IgG, the IgM quotient Q IgM and the albumin quotient Q Alb are calculated according to the following formulas:

- (4) IgG-quotient $Q_{\text{IgG}} = \text{IgG conc. CSF} / \text{IgG IgG conc. serum}$
- (5) IgM-quotient $Q_{\text{IgM}} = \text{IgM conc. CSF} / \text{IgM IgG conc. serum}$
- (6) albumin quotient $Q_{\text{Alb}} = \text{albumin conc. CSF} / \text{albumin conc. serum}$

5.2.4. Calculation of the limiting values QLim (IgG) and QLim (IgM)

- (7) $Q_{\text{Lim}}(\text{IgG}) = 0,93 \times \sqrt{Q_{\text{Alb}}^2 + (6 \times 10^{-6})} - 1,7 \times 10^{-3}$
- (8) $Q_{\text{Lim}}(\text{IgM}) = 0,67 \times \sqrt{Q_{\text{Alb}}^2 + (120 \times 10^{-6})} - 7,1 \times 10^{-3}$

5.2.5. Calculation of the antibody index (AI)

The values obtained for QIgG and QLim (IgG) are compared and AI (IgG) is calculated using formula (9) or (10). Correspondingly, QIgM and QLim (IgM) are compared and AI (IgM) is calculated using formula (11) or (12).

if $QIgG < QLim (IgG)$

$$(9) \quad \mathbf{AI (IgG)} = Qspec (IgG) / QIgG$$

If $QIgG > QLim (IgG)$

$$(10) \quad \mathbf{AI (IgG)} = Qspec (IgG) / QLim (IgG)$$

if $QIgM < QLim (IgM)$

$$(11) \quad \mathbf{AI (IgM)} = Qspec (IgM) / QIgM$$

If $QIgM > QLim (IgM)$

$$(12) \quad \mathbf{AI (IgM)} = Qspec (IgM) / QLim (IgM)$$

5.2.6 Interpretation of AI values

$0,6 < AI (IgG) \text{ resp. } (IgM)$	normal
$1,3 < AI < 1,5 (IgG) \text{ resp. } (IgM)$	borderline
$AI (IgG) \text{ resp. } (IgM) \geq 1,5$	pathological

Antibody index values < 0.6 indicate an analytical fault (e.g. clinical chemistry). It is recommended to check whether serum and CSF are taken from the same day and if sample-mix-up can be excluded. A repetitive determination of the values of clinical chemistry and/ or of the Borrelia specific antibodies might be advisable.

Antibody index values $0.6 < AI < 1.3$ are in the normal range. An intrathecal synthesis of antibodies against Borrelia is unlikely. The detected Borrelia-antibodies in the CSF originate from the serum and have passively passed into the CSF space through the blood-liquor barrier.

Antibody index values between 1.3 and 1.5 are borderline cases. The serum-liquor pair should be retested. Or if clinical doubt exists, it is recommended to test another serum-CSF pair in the course of the infection.

Antibody index values ≥ 1.5 are to be classified as pathological. They indicate intrathecal synthesis of Borrelia-antibodies as found in patients suffering from neuroborreliosis.

6. SPECIFIC PERFORMANCE CHARACTERISTICS

6.1. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

It is $> 90 \%$.

6.2. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

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NovaTec Immundiagnostica GmbH