Leishmania Canine

IgG - ELISA

Enzyme immunoassay for the qualitative detection of IgG-class antibodies against Leishmania in canine serum
Enzymimmunoassay zum qualitativen Nachweis von IgG-Antikörpern gegen Leishmania in Hundeserum

Only for in-vitro diagnostic use

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ENGLISH

1. INTRODUCTION

Leishmania are protozoa belonging to the family Trypanosomatidae. They are dimorphic obligate intracellular parasites which are transmitted as flagellated forms (promastigotes) by the bite of several species of phlebotomine sandflies. Once inoculated into the skin of the mammalian host, the promastigotes are engulfed by macrophages, in which they transform to amastigotes, an aflagellated form. Inside the macrophages the amastigotes divide and spread to different organs, especially the bone marrow, lymph nodes, skin, spleen, liver, and kidneys.

Leishmania infantum infected dogs constitute the main reservoir of the parasite. Canine visceral leishmaniasis is widespread in countries of the Mediterranean basin, Asia and Latin America. In addition, imported cases of infection are occasionally reported in Northern Europe and North America.

Exposure to bites of infected sandflies can result in an asymptomatic period that may progress to symptomatic infection or remain cryptic. However, both symptomatic and asymptomatic animals harbor the parasite and are infective to sandfly vectors.

Clinical symptoms may appear months to several years after infection. The clinical picture varies widely; symptoms of leishmaniasis include skin lesions (dry exfoliative dermatitis, ulcerations and alopecia), weight loss, local or generalized lymphadenopathy, ocular and periocular lesions, renal failure, epistaxis, lameness, and anemia. Occasionally, some dogs develop chronic diarrhea or liver failure. Polyclonal hyperproteinemia is observed in most dogs. Some animals show leukopenia, whereas others develop leukocytosis. In animals with renal lesions, increased plasma urea and creatinine, proteinuria, and hematuria are common.

2. INTENDED USE

The NovaTec Leishmania Canine IgG-ELISA is intended for the qualitative determination of IgG-class antibodies against Leishmania in canine serum.

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of IgG-class antibodies against Leishmania is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiter strip wells are precoated with Leishmania antigens to bind corresponding antibodies of the specimen. After washing the wells to remove all unbound sample material horseradish peroxidase (HRP) labelled anti-dog IgG conjugate is added. This conjugate binds to the captured Leishmania-specific antibodies. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine- (TMB) substrate which gives a blue reaction product. The intensity of this product is proportional to the amount of Leishmania-specific IgG antibodies in the specimen. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Leishmania Microtiter Strips (IgG):** 12 breakapart 8-well snap-off strips coated with Leishmania antigen; in resealable aluminium foil.
- **IgG Sample Diluent***:** 1 bottle containing 100 ml of buffer for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulfuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Solution (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated buffer (pH 7.2 ± 0.2) for washing the wells; white cap.
- **Leishmania anti-dog IgG Conjugate**:** 1 bottle containing 20 ml of peroxidase labelled anti-dog IgG, coloured red, ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB); ready to use; yellow cap.
- **Leishmania IgG Positive Control***:** 1 bottle containing 2 ml; coloured yellow; ready to use; red cap.
- **Leishmania IgG Negative Control***:** 1 bottle containing 2 ml; coloured yellow; ready to use; blue cap.
- **Leishmania IgG Cut-off Control***:** 1 bottle containing 3 ml, coloured yellow, ready to use, green cap.

* contains 0.01 % Kathon after dilution

** contains 0.2 % Bronidox L

*** contains 0.1 % Kathon

4.2. Materials supplied

- 1 Strip holder
- 2 Cover foils
- 1 Test protocol
- 1 distribution and identification plan

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at +2...+8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and controls to room temperature (+20...+25 °C) before starting the test run!

6.1. Microtiter Strips

The ready to use breakapart snap-off strips are coated with Leishmania antigen. Store at +2...+8 °C. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at +2...+8 °C; stability until the expiry date.*

6.2. Leishmania anti-dog IgG Conjugate

The bottle contains 20 ml of a solution with anti-dog IgG conjugated to horseradish peroxidase, buffer, stabilizers, preservatives and an inert red dye. The solution is ready to use. Store at +2...+8 °C. *After first opening stability until the expiry date when stored at +2...+8 °C.*

6.3. Controls

The bottles contain a ready to use control solution. They have to be stored at +2...+8 °C. *After first opening stability until the expiry date when stored at +2...+8 °C.*

6.4. IgG Sample Diluent

The bottle contains 100 ml phosphate buffer, stabilizers, preservatives and an inert yellow dye. It is used for the dilution of the patient specimen. This ready to use solution has to be stored at +2...+8 °C. *After first opening stability until the expiry date when stored at +2...+8 °C.*

6.5. Washing Solution (20x conc.)

The bottle contains 50 ml of a concentrated buffer, detergents and preservatives. Dilute washing solution 1+19; e.g. 10 ml washing solution + 190 ml fresh and germ free redistilled water. The diluted buffer is stable for 5 days if stored at room temperature. *Crystals in the solution disappear by warming up to +37 °C in a water bath. After first opening the concentrate is stable until the expiry date.*

6.6. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at +2...+8 °C, away from the light. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away. After first opening stability until the expiry date when stored at +2...+8 °C.*

6.7. Stop Solution

The bottle contains 15 ml 0.2 M sulphuric acid solution (R 36/38, S 26). This ready to use solution has to be stored at +2...+8 °C. *After first opening stability until the expiry date.*

7. SPECIMEN COLLECTION AND PREPARATION

Use dog serum samples with this assay. If the assay is performed within 5 days after sample collection, the specimens should be kept at +2...+8 °C; otherwise they should be aliquoted and stored deep-frozen (-20 to -70 °C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Diluent. Dispense 10 µl sample and 1 ml IgG Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex. *Controls are ready to use and must not be diluted.*

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described. If performing the test on ELISA automatic systems we recommend to increase the washing steps from three to five and the volume of washing solution from 300 µl to 350 µl to avoid washing effects. Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

1 well	(e.g. A1)	for the substrate blank,
1 well	(e.g. B1)	for the Negative Control,
2 wells	(e.g. C1+D1)	for the Cut-off Control and
1 well	(e.g. E1)	for the Positive Control.

It is recommended to determine controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each control and sample.

Adjust the incubator to $+37 \pm 1$ °C.

1. Dispense 100 µl controls and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 60 min \pm 5 min at $+37 \pm 1$ °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of Washing Solution. Avoid overflows from the reaction wells. The soak time between each wash cycle should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µl Leishmania anti-dog IgG Conjugate into all wells except for the blank well (e.g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature. Do not expose to direct sunlight.**
7. Repeat step 4.
8. Dispense 100 µl TMB Substrate Solution into all wells
9. **Incubate for exactly 15 min at room temperature in the dark.**
10. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution.
Any blue colour developed during the incubation turns into yellow.
11. Measure the absorbance of the specimen at 450/620 nm within 30 min after addition of the Stop Solution.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each control and patient sample in the distribution and identification plan.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate blank** in A1: Absorbance value **lower than 0.100**
- **Negative Control** in B1: Absorbance value **lower than 0.200**
- **Cut-off Control** in C1 and D1: Absorbance value between **0.250 and 0.900**
- **Positive Control** in E1: Absorbance value higher than Cut-off

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off value is calculated as the mean value of the determined Cut-off controls in C1 and D1.

Example: Absorbance value Cut-off Control 0.39 + absorbance value Cut-off Control 0.37 = $0.76 / 2 = 0.38$

Cut-off = 0.38

9.3. Interpretation of Results

Samples are considered **POSITIVE** if the absorbance value is higher than 20 % over the cut-off.

Samples with an absorbance value of 20 % above or below the cut-off should not be considered as clearly positive or negative
→ **grey zone**

It is recommended to repeat the test again 2 - 4 weeks later with a fresh sample. If results in the second test are again in the grey zone the sample has to be considered **NEGATIVE**.

Samples are considered **NEGATIVE** if the absorbance value is lower than 20 % below the cut-off.

9.3.1. Results in NovaTec Units

$$\frac{\text{Patient (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec-Units} = \text{NTU}]$$

Example:
$$\frac{1.786 \times 10}{0.43} = 42 \text{ NTU (NovaTec Units)}$$

Cut-off:	10	NTU
Grey zone:	8-12	NTU
Negative:	<8	NTU
Positive:	>12	NTU

10. SPECIFIC PERFORMANCE CHARACTERISTICS

10.1. Precision

Interassay	n	Mean	Cv (%)
Pos. Serum	4	1.54	7.1

Intraassay	n	Mean	Cv (%)
Pos. Serum	11	1.6	6.1

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.
It is > 95 %.

10.3 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.
It is 95 %.

Note: The results refer to the groups of samples investigated; these are not guaranteed specifications.
--

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values. Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after the expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The NovaTec ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

WARNING:	In the used concentration Bronidox L has hardly any toxicological risk upon contact with skin and mucous membranes!
WARNING:	Sulphuric acid irritates eyes and skin. Keep out of the reach of children. Upon contact with the eyes, rinse thoroughly with water and consult a doctor!

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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
Leishmania Canine IgG-ELISA

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all specimens and controls on the result sheet supplied in the kit.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank (e.g. A1)	Negative Control	Positive Control	Cut-off Control	Sample (diluted 1+100)
Negative Control	-	100µl	-	-	-
Positive Control	-	-	100µl	-	-
Cut-off Control	-	-	-	100µl	-
Sample (diluted 1+100)	-	-	-	-	100µl
Cover wells with foil supplied in the kit Incubate for 1 h at 37°C Wash each well three times with 300µl of washing solution					
Conjugate	-	100µl	100µl	100µl	100µl
Cover wells with foil supplied in the kit Incubate for 30 min at room temperature Wash each well three times with 300µl of washing solution					
TMB Substrate	100µl	100µl	100µl	100µl	100µl
Incubate for exactly 15 min at room temperature in the dark					
Stop Solution	100µl	100µl	100µl	100µl	100µl
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

NovaTec Immundiagnostica GmbH

Technologie & Waldpark

Waldstr. 23 A6
 D-63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760
 Email: info@NovaTec-ID.com
 Internet: www.NovaTec-ID.com

Fax: +49 (0) 6074-487629

LEIG0310V-canine-dt,engl-26082008-EH