Borrelia burgdorferi

Lyme Blot Canine IgG

Immunoblot for the detection of IgG-class antibodies against *Borrelia burgdorferi* sensu lato in canine serum

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ENGLISH

1. INTRODUCTION

Lyme disease is a tick-borne multistage disease which is caused by spirochetal bacteria of the *Borrelia burgdorferi* sensu lato complex. In Europe *Borrelia afzelii* and *Borrelia garinii* are the most frequently isolated pathogenic genomic species in man and dog. *Borrelia burgdorferi* sensu stricto is known to play an inferior role in Europe.

Borrelia infections in the dog are quite common but often occur asymptotically. Clinical symptoms are strongly dependent on the particular *Borrelia* species. Clinical manifestations of canine Lyme disease, including fever, anorexia, lethargy, lymphadenopathia and alternating lameness due to polyarthritis are mainly attributed to the species *B.b. sensu stricto*. The characteristic rash or the circular area of redness around the bite (erythema chronicum migrans) which is seen in man is often absent in the dog. Sometimes a small reddened patch occurs which may be overlooked due to hair coat or dark pigmentation and frequently disappears within the first week.

Today, the laboratory workup of Lyme disease relies on a good anamnesis, an extensive medical checkup as well as purpose-oriented laboratory diagnostics. The standard approach involves screening the patient's sample using a technology such as ELISA or IFA. Due to false-positive and false-negative test results a positive result should be verified by a confirmatory test like the NovaTec *Borrelia burgdorferi* Lyme Blot Canine IgG.

The immune reaction to a *Borrelia* infection is complex as the delay in the development of detectable antibodies may vary in different individuals. Antibody production can be affected by various criteria like *Borrelia* genotype, rate of infection, age and immune status of the dog, medical treatment with glucocorticoids etc. Cross-reactivity of some *Borrelia* proteins with antigens from other bacteria is also well-known. The Western Blot technology is extremely useful in dissecting the immune response to *Borrelia* infections which develops gradually over a period of weeks to years and which involves the appearance of IgM and IgG antibodies directed against a number of *Borrelia*-associated proteins.

2. INTENDED USE

The NovaTec *Borrelia burgdorferi* Lyme Blot Canine IgG is a qualitative assay for the detection of *Borrelia burgdorferi* sensu lato specific IgG antibodies in canine serum.

3. PRINCIPLE OF THE ASSAY

Proteins derived from *Borrelia garinii* were separated electrophoretically according to their molecular weight using sodium dodecyl sulfate polyacrylamide gel electrophoresis (SDS-PAGE). Following electrophoresis, the separated protein bands were transferred to a PVDF membrane. The antigen bearing membrane was blocked and cut into ready-to-use strips.

Diluted dog sera are incubated with the Antigen Strips. *Borrelia*-specific antibodies, if present, will bind to their target antigens. After washing, anti-dog IgG conjugated with horseradish peroxidase is added. At the end of a second incubation, unbound conjugate is removed by washing and aspiration. The bound conjugate is visualized by the addition of a chromogenic substrate. Strips are dried and analysed. Using the kit-specific Template supplied with the kit, the position of the stained bands can be correlated with defined *Borrelia* antigens.

4. MATERIALS

4.1. Reagents supplied

- **Borrelia garinii Antigen Strips (IgG):** 1 tube containing 10 consecutively numbered western blot strips coated with *Borrelia garinii* antigens.
- **Template:** 1 predeveloped control strip, kit-specific.
- **IgG Sample Diluent***:** 1 bottle containing 50 ml of buffer for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- **Washing Solution (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated buffer, pH 7.2 ± 0.2 for washing the Antigen Strips; white cap.
- **Lyme anti-dog IgG Conjugate**:** 1 bottle containing 12 ml of peroxidase labelled anti-dog IgG; coloured red, ready to use; black cap.
- **TMB Substrate-Solution Membrane:** 1 bottle containing 12 ml 3,3',5,5'-tetramethylbenzidine (TMB); ready to use; blue cap.

* contains 0.1 % Bronidox L after dilution

** contains 0.2 % Bronidox L

*** contains 0.1 % Kathon

4.2. Materials supplied

- 1 test protocol
- 1 result sheet
- 2 incubation trays

4.3. Materials and Equipment needed

- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Shaking platform
- Micropipettes, 10 and 1000 µl
- Plastic tweezers for handling Antigen Strips
- Vacuum apparatus
- Timer
- Filter paper

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at +2...+8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and controls to room temperature (+20...+25 °C) before starting the test run!

6.1. Borrelia garinii Antigen Strips

The ready to use western blot strips are coated with Borrelia garinii antigens. Store at +2...+8 °C. *After first opening stability until the expiry date when stored at +2...+8 °C. Do not interchange strips from different kits!*

6.2. Lyme anti-dog IgG Conjugate

The bottle contains 12 ml of a solution with anti-dog IgG conjugated to horseradish peroxidase, buffer, stabilizers, preservatives and an inert red dye. The solution is ready to use. Store at +2...+8 °C. *After first opening stability until the expiry date when stored at +2...+8 °C.*

6.3. IgG Sample Diluent

The bottle contains 50 ml phosphate buffer, stabilizers, preservatives and an inert yellow dye. It is used for the dilution of the patient specimen. This ready to use solution has to be stored at +2...+8 °C. *After first opening stability until the expiry date when stored at +2...+8 °C.*

6.4. Washing Solution (20x conc.)

The bottle contains 50 ml of a concentrated buffer, detergents and preservatives. Dilute washing solution 1+19; e.g. 10 ml Washing Solution + 190 ml fresh and germ free redistilled water. The diluted buffer is stable for 5 days if stored at room temperature. *Crystals in the solution disappear by warming up to +37 °C in a water bath. After first opening the concentrate is stable until the expiry date when stored at +2...+8 °C.*

6.5. TMB Substrate-Solution Membrane

The bottle contains 12 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at +2...+8 °C, away from the light. *The solution should be colourless to pale yellow. If the substrate turns into blue, it may have become contaminated and should be thrown away. After first opening stability until the expiry date when stored at +2...+8 °C.*

7. SPECIMEN COLLECTION AND PREPARATION

Use dog serum samples with this assay. If the assay is performed within 5 days after sample collection, the specimens should be kept at +2...+8 °C; otherwise they should be aliquoted and stored deep-frozen (-20...-70 °C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing!*

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Diluent. Dispense 12 µl sample and 1.2 ml IgG Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described. Perform all assay steps in the correct order and without any appreciable delays between the steps. A clean, disposable tip should be used for dispensing each sample.

1. Using tweezers, carefully remove the number of required strips from the tube and place them in the channels of the incubation tray.
2. Transfer 1.2 ml of diluted dog serum to the appropriate channel of the incubation tray. A new pipette tip must be used with each separate sample. Check that all Antigen Strips are **completely immersed** and, if needed, gently shake the tray or push delicately the strips into the solution with a clean pipette tip. The **numbered side of the strip must face up**; antigens are bound to this side of the membrane.

3. Incubate the strips at room temperature for **60 min** on a shaking platform.
4. Wash procedure
 - Carefully aspirate the contents of each channel, add 1 ml of Washing Solution.
 - Place the incubation tray on the shaking platform for 5 min at room temperature.
 - Aspirate the contents of each channel and repeat the wash procedure two more times.
5. Add **1 ml** of anti-IgG Conjugate to the appropriate channels and incubate at room temperature for **30 min** on the shaking platform.
6. Repeat step 4.
7. Add **1 ml** of chromogenic TMB Substrate-Solution Membrane to the appropriate channels and incubate for **10 min** at room temperature and under observation on the shaking platform.
CAUTION: DO NOT OVERDEVELOP. Some sera show excessive background staining of the strips. To avoid background, stop the reaction earlier by washing with deionised water.
8. Stop the reaction by aspirating the contents of each channel and rinse each strip with **1 ml** of deionised water. Place back the incubation tray on the shaking platform for 5 min. Repeat the 1 ml - 5 min rinse.
9. Remove the strips from the incubation tray and place them with the number face up on filter paper to dry (approximately 30 min at room temperature).

IMPORTANT: to prevent fading, the Antigen Strips should be protected from exposure to light.

9. RESULTS

9.1. Interpretation of Results

Once dried, attach the Antigen Strips to the result sheet using clear tape.

- The *Borrelia burgdorferi* Lyme Blot Canine IgG utilizes a kit-specific Template consisting of a developed western blot strip cut from the membrane used to prepare the Antigen Strips. This strip has been exposed to a positive control serum in order to exhibit bands representative of a Lyme infection. The corresponding protein bands are listed on the right side of the Template.
- Identification of the reactive bands is based upon comparison with the exposed bands on the Template. The horizontal line at the bottom of the strip must be aligned with the index line near the bottom of the control strip.
- For band scoring, determine the reactivity of the following bands:
p100 (93 kDa), **p58** (58 kDa), **flagellin** (41 kDa), **BmpA** (39 kDa), **OspA** (32.5 kDa), **OspC** (22-23 kDa) and **p18** (18 kDa)

9.2. Band Specificity

Band	Antigen Characterization and their diagnostic Relevance
p100 (93 kDa)	chromosomal protein highly specific of the <i>Borrelia</i> genus; the antibody response is mostly IgG and appears in the course of chronic infections
75 kDa	specific marker, probably corresponding to the bacterial DNA
62-72 kDa	heat shock proteins; these bands are not specifically associated with a <i>Borrelia</i> infection and are found in several bacterial infections; no diagnostic value
66 kDa	species-dependent but specific marker for <i>Borrelia</i> infections
60 kDa	common bacterial antigen; non-specific; no diagnostic value
58 kDa	species-dependent but specific marker for <i>Borrelia</i> infections => infections with high <i>Borrelia</i> load (often in combination with p18) => after vaccinations (relatively specific vaccination marker)
flagellin (41 kDa)	Flagellin protein => cross-reaction with other spirochetes are common => all infection-stages
BmpA (39 kDa)	highly specific Borreliac membrane protein A; detectable at all stages of infection
37 kDa	this protein is considered to be an early marker but its specificity has not been established

OspA (32.5 kDa)	Outer surface protein A; highly specific of <i>Borrelia garinii</i> => after infection: in later stages, low intensity, rarely (in case of multiple antigen contacts stronger) => after vaccination: early, very intensive and long-lasting
OspC (22-23 kDa)	Outer surface protein C; highly specific of <i>Borrelia garinii</i> => after infection: often in early stages (important IgM-marker), but also in late stages => after vaccination: occasionally (early and short appearance)
18 kDa	Highly specific ; excellent marker for late infections (particularly in strong positive sera)

Dog sera may exhibit other bands than those mentioned above. Such bands should not be considered when interpreting the results.

9.3. Evaluation of Antigen Bands

- Before analyzing the bands, two questions have to be addressed:
 1. Is the animal frequently exposed to tick bites? If this is the case, the blot will automatically show several specific bands even in the absence of an ongoing infection. Other serological tests will probably be positive and a complete analysis of the Western blot bands will be needed to differentiate between a past exposure and an ongoing infection. In contrast, dogs which are not exposed to tick bites will show immunoblots, when negative, with very few bands (most frequent bands: 60 and 41 kDa).
 2. Is the clinical presentation suggestive of an ongoing borreliosis? In the absence of fever and/or lameness, distinguishing between a past exposure and an ongoing infection may be extremely difficult.
- Overall presentation of the Antigen Strip: a positive result corresponds to a strip with **several neat bands**, including some very intense bands. Positive blots are very similar to the control strip on the kit Template.

IMPORTANT: only bands with the same or stronger color intensities than the according band on the kit-specific control strip (Template) are considered for evaluation.

Point evaluation of *Borrelia* antigens in the *Borrelia burgdorferi* Lyme Blot Canine IgG

Borrelia Antigen	Molecular Weight [kDa]	Points IgG
p100	93	5
p58	58	3
flagellin	41	1
BmpA	39	5
OspA	32.5	3
OspC	22-23	5
p18	18	5

The following important points form the basis of the final sum up:

1. No blot result based on a single immunodominant band is to be considered as positive
2. If only one of the following bands: p100, p39, OspC or p18 is reactive, a second specimen should be tested after 2-3 weeks using the *Borrelia burgdorferi* Lyme Blot Canine IgG
3. As soon as at least 2 bands (p100, p39, OspC or p18) are considered to be positive, the likelihood of a natural *Borrelia* infection is > 95 %

There are no classical “vaccination” or “infection bands”, but there are bands reacting stronger after vaccination than after infection (e.g. OspA).

Sum up the points for each patient specimen and interpret as follows:

Sum of points	Evaluation IgG
< 5	negative
5-7	equivocal
≥ 8	positive

10. LIMITATIONS OF THE PROCEDURE

The absence of IgG antibodies reactive with *Borrelia* antigens on the *Borrelia burgdorferi* Lyme Blot Canine IgG cannot exclude a *Borrelia* infection in all cases. In the early phase of infection, detectable amounts of antibodies may not or not yet be present. In addition, antibiotic therapy in early stages of infection can suppress the formation of antibodies.

The presence of antibodies specific for *Borrelia* antigens does not necessarily indicate an acute infection since high IgG titers of a previous infection can possibly persist for several years.

Serological results should not be used as the sole criterion on which a diagnosis is established. They are an aid to diagnosis only, and their significance must be interpreted in relationship to the clinical presentation of the animal. In the case of unclear or equivocal findings, repeated testing extending over several weeks is recommended.

11. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons.
- Only for in-vitro diagnostic use.
- All materials should be regarded and handled as potentially infectious.
- Do not interchange reagents from different production lots or Antigen Strips from different kits.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after the expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check reagent vials for microbial contamination prior to further use.
- Always use tweezers to handle Antigen Strips.
- During incubations and washing, the Antigen Strips must remain completely covered with fluid and the numbered side of the strips must face up.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The NovaTec Blot is only designed for qualified personnel who are familiar with good laboratory practice.

WARNING:	In the used concentration Bronidox L has hardly any toxicological risk upon contact with skin and mucous membranes.
-----------------	---

11.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Hovius J.W.R., Hovius K.E., Oei A., Houwers D.J., van Dam A.P. (2000) Antibodies against Specific Proteins of and Immobilizing Activity against Three Strains of *Borrelia burgdorferi* Sensu Lato Can Be Found in Symptomatic but Not in Infected Asymptomatic Dogs. J Clin Microbiol 38, 2611-2621

Guerra M.A., Walker E.D., Kitron U. (2000) Quantitative Approach for the Serodiagnosis of Canine Lyme Disease by the Immunoblot Procedure. J Clin Microbiol 38, 2628-2632

Symbols Key / Symbolschlüssel	
	Manufactured by / Hergestellt von
IVD	In Vitro Diagnostic Medical Device / In Vitro Diagnosticum
LOT	Lot Number / Chargenbezeichnung
	Expiration Date / Verfallsdatum
	Storage Temperature / Lagertemperatur
[REF]	Catalogue Number / Katalog Nummer
	Consult Instructions for Use / Gebrauchsanweisung beachten
STRIP	Antigen-Strips / Antigen-Streifen
TEMP	Template / Schablone
CONJ	Conjugate / Konjugat
DIL G	Sample diluent buffer IgG / IgG-Probenverdünnungspuffer
SUB TMB	TMB Substrate-Solution Membrane / TMB-Substrat-Lösung Membran
WASH BUF 20x	Washing solution 20x concentrated / Waschlösung 20x konzentriert
	Contains sufficient for "n" tests / Ausreichend für "n" Tests

SCHEME OF THE ASSAY

Borrelia burgdorferi Lyme Blot Canine IgG

Test Preparation

Prepare samples as described.
Note sample and strip numbers as well as lot number.
Put the required number of Antigen Strips into the wells of the incubation tray.
Prepare the required volume of Washing Solution.

Assay Procedure

pipette 1.2 ml diluted sample
↓
incubate for 60 min at room temperature while shaking gently
↓
wash three times for 5 min with 1 ml Washing Solution while shaking
↓
add 1 ml Conjugate
↓
incubate for 30 min at room temperature while shaking gently
↓
wash three times for 5 min with 1 ml Washing Solution while shaking
↓
add 1 ml Substrate-Solution
↓
incubate for 10 min at room temperature while shaking gently
↓
wash two times for 5 min with 1 ml deionised water
dry and read

NovaTec Immundiagnostica GmbH

Technologie & Waldpark

Waldstr. 23 A6
D-63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760
Email: info@NovaTec-ID.com
Internet: www.NovaTec-ID.com

Fax: +49 (0) 6074-487629

LYGV0110engl,dt080609-EH