## NovaLisa™

# Parainfluenza Virus 1, 2, 3

## IgA – ELISA

Enzyme immunoassay for the qualitative determination of IgA-class antibodies against Parainfluenza Virus (Type 1, 2 and 3) in human serum or plasma

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For further languages please contact our authorized distributors.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografia	Page / Seite / Page / Pagina / Página	26
Symbols Key / Symbolschlüssel / Explication des symboles / Legenda / Símbolos	Page / Seite / Page / Pagina / Página	27
Summary of Test Procedure/ Kurzanleitung Testdurchführung/ Résumé de la procédure de test/ Schema della procedura/ Resumen de la técnica	Page / Seite / Page/ / Pagina / Página	28

## 1. INTRODUCTION

Parainfluenza viruses are important viral pathogens causing upper and lower respiratory infections in adults and children. The viruses belong to the family of Paramyxoviridae. They are enveloped single-stranded RNA viruses with spherical or pleomorphic shape. Parainfluenza viruses are relatively large viruses of about 150-300 nm in diameter. On the basis of antigenic differences they are divided into four subtypes, of which type 4 is divided into two subtypes. Parainfluenza virus 1 and 3 belong to the paramyxovirus genus. Parainfluenza virus 2, 4a and 4b belong to the rubella virus genus along with mumps. The virus is ubiquitous; infections occur as epidemics as well as sporadically. Parainfluenza viruses are sensitive to detergents and heat but can remain viable on surfaces for up to 10 hours. Transmission occurs via droplets, aerosols and fomites (viruses survive on surfaces). Parainfluenza virus infections are primarily childhood diseases, the highest age-specific attack rates for croup occur in children below the age of 3 years. Type 3 infections occur earliest and most frequently, so that ~50% of children are infected during the first year of life and almost all by 6 years, as determined by seroepidemiological studies. Antibodies against types 1 and 2 are acquired less rapidly but 80% of children have antibodies by 10 years of age. Type 4 viruses induce few clinical illnesses but infections are common, 70 - 80% of children have antibodies by 10 years of age. In all age groups, the incubation period appears to be five to six days. Croup or laryngotracheobronchitis is the commonest clinical manifestation of infection. Parainfluenza viruses are found uncommonly associated with other respiratory tract infections in children such as tracheobronchitis, bronchiolitis, and bronchopneumonia. Occasionally, a mild non-specific illness is seen after parainfluenza virus infection. In adults, the virus is usually limited to causing inflammation in the upper parts of the respiratory tract.

Species	Disease	Symptoms	Mechanism of Infection
Parainfluenza Virus (Type 1, 2 and 3)	The majority of parainfluenza virus infections result in a non-specific upper respiratory tract infection, even in children. Parainfluenza virus infections in adults are relatively uncommon, and symptoms are usually less severe in adults than in children.	Milder, influenza-like illness; fever, headache, cough, sore throat; bronchitis.  Complications: bronchopneumonia, croup, bacterial superinfection	Infections are spread via: droplets, aerosols (upon contact) fomites (indirect contact transmission)

Infections can be confirmed in two ways:

- by isolation and identification of the virus in cell culture or by direct detection of the virus in respiratory secretions (usually, collected within one week of onset of symptoms) using immunofluorescence, enzyme immunoassay, or polymerase chain reaction assay, and
- by demonstration of a significant rise in specific IgA antibodies between appropriately collected paired serum specimens or specific IgM antibodies in a single serum specimen

## 2. INTENDED USE

The NovaTec Parainfluenza Virus Type 1, 2, 3 IgA-ELISA is intended for the qualitative determination of IgA class antibodies against Parainfluenza Virus IgA in human serum or plasma (citrate).

## 3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of IgA-class antibodies against Parainfluenza Viruses (Type 1, 2 and 3) is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microwells are precoated with Parainfluenza Virus antigens to bind corresponding antibodies of the specimen. After washing the wells to remove all unbound sample material horseradish peroxidase (HRP) labelled anti-human IgA conjugate is added. This conjugate binds to the captured Parainfluenza Virus-specific antibodies. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is proportional to the amount of Parainfluenza Virus-specific IgA antibodies in the specimen. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450 nm is read using an ELISA microwell plate reader.

## 4. MATERIALS

### 4.1. Reagents supplied

- Parainfluenza Virus Coated Wells (IgA):** 12 breakapart 8-well snap-off strips coated with Parainfluenza Virus antigen; in resealable aluminium foil.
- IgA Sample Diluent \*\*\*:** 1 bottle containing 100 ml of buffer for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- Washing Solution (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated buffer (pH 7.2 ± 0.2) for washing the wells; white cap.
- Parainfluenza Virus anti-IgA Conjugate\*\*\*:** 1 bottle containing 20 ml of peroxidase labelled rabbit antibody to human IgA; coloured violet, ready to use; black cap.

- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB); ready to use; yellow cap.
  - **Parainfluenza Virus IgA Positive Control\*\*\*:** 1 bottle containing 2 ml; coloured yellow; ready to use; red cap.
  - **Parainfluenza Virus IgA Cut-off Control\*\*\*:** 1 bottle containing 3 ml; coloured yellow; ready to use; green cap.
  - **Parainfluenza Virus IgA Negative Control\*\*\*:** 1 bottle containing 2 ml; coloured yellow; ready to use; blue cap.
- \* contains 0.1 % Bronidox L after dilution  
 \*\* contains 0.2 % Bronidox L  
 \*\*\* contains 0.1 % Kathon

## 4.2. Materials supplied

- 1 Strip holder
- 1 Cover foil
- 1 Test protocol
- 1 distribution and identification plan

## 4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Timer

## 5. STABILITY AND STORAGE

---

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

## 6. REAGENT PREPARATION

---

*It is very important to bring all reagents, samples and controls to room temperature (20...25°C) before starting the test run!*

### 6.1. Coated snap-off Strips

The ready to use breakapart strips are coated with Parainfluenza Virus antigen. Store at 2...8°C. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

### 6.2. Parainfluenza Virus anti-IgA Conjugate

The bottle contains 20 ml of a solution with anti-human-IgA horseradish peroxidase, buffer, stabilizers, preservatives and an inert violet dye. The solution is ready to use. Store at 2...8°C. *After first opening stability until expiry date when stored at 2...8°C.*

### 6.3. Controls

The bottles labelled with Positive, Cut-off and Negative Control contain a ready to use control solution. It contains 0.1% Kathon and has to be stored at 2...8°C. *After first opening stability until expiry date when stored at 2...8°C.*

### 6.4. IgA Sample Diluent

The bottle contains 100 ml phosphate buffer, stabilizers, preservatives and an inert yellow dye. It is used for the dilution of the patient specimen. This ready to use solution has to be stored at 2...8°C. *After first opening stability until expiry date when stored at 2...8°C.*

### 6.5. Washing Solution (20xconc.)

The bottle contains 50 ml of a concentrated buffer, detergents and preservatives. Dilute washing solution 1+19; e.g. 10 ml washing solution + 190 ml fresh and germ free redistilled water. The diluted buffer will keep for 5 days if stored at room temperature. *Crystals in the solution disappear by warming up to 37 °C in a water bath. After first opening the concentrate is stable until the expiry date*

### 6.6. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C, away from light. *The solution should be colourless or have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be discharged. After first opening stability until expiry date when stored at 2...8°C.*

### 6.7. Stop Solution

The bottle contains 15 ml 0.2 M sulphuric acid solution (R 36/38, S 26), ready to use, store at 2...8°C. *After first opening stability until expiry date.*

## 7. SPECIMEN COLLECTION AND PREPARATION

---

Use human serum or plasma (citrate) samples with this assay. If the assay is performed within 5 days after sample collection, the specimen should be kept at 2...8°C; otherwise they should be aliquoted and stored deep-frozen (-20 to -70°C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*  
Heat inactivation of samples is not recommended.

### 7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgA Sample Diluent. Dispense 10µl sample and 1 ml IgA Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

## 8. ASSAY PROCEDURE

---

### 8.1. Test Preparation

Please read the test protocol carefully before performing the assay. Result reliability depends on strict adherence to the test protocol as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend to increase the washing steps from three to five and the volume of washing solution from 300µl to 350µl to avoid washing effects. Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

1 well	(e.g. A1)	for the substrate blank,
1 well	(e.g. B1)	for the negative control,
2 wells	(e.g. C1+D1)	for the cut-off control and
1 well	(e.g. E1)	for the positive control.

*It is recommended to determine controls and patient samples in duplicate.*

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each control and sample.

Adjust the incubator to 37° ± 1°C.

1. Dispense 100µl controls and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37±1°C.**
4. When incubation has been completed, remove the foil, aspirate the content off the wells and wash each well three times with 300µl of Washing Solution. Avoid overflows from the reaction wells. The soak time between each wash cycle should be >5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!  
*Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.*
5. Dispense 100µl Parainfluenza Virus anti IgA Conjugate into all wells except for the blank well (e.g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature. Do not expose to direct sunlight.**
7. Repeat step 4.
8. Dispense 100µl TMB Substrate Solution into all wells
9. **Incubate for exactly 15 min at room temperature in the dark.**
10. Dispense 100µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution.  
*Any blue colour developed during the incubation turns into yellow.*

*Note: Highly positive patient samples can cause dark precipitates of the chromogen! These precipitates have an influence when reading the optical density. Predilution of the sample with physiological sodium chloride solution, for example 1+1, is recommended. Then dilute the sample 1+100 with dilution buffer and multiply the results in NTU by 2.*

11. Measure the absorbance of the specimen at 450/620 nm within 30 min after addition of the Stop Solution.

### 8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

*If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!*

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each control and patient sample in the distribution and identification plan.

*Dual wavelength reading using 620 nm as reference wavelength is recommended.*

Where applicable calculate the **mean absorbance values** of all duplicates.

## 9. RESULTS

---

### 9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate blank** in A1: Absorbance value < **0.100**.
- **Negative control** in B1: Absorbance value < **0.200 and < cut-off**
- **Cut-off control** in C1 and D1: Absorbance value **0.150 – 1.30**.
- **Positive control** in E1: Absorbance value > **cut-off**.

If these criteria are not met, the test is not valid and must be repeated.

### 9.2. Calculation of Results

**The cut-off is the mean absorbance value of the Cut-off control determinations.**

*Example: Absorbance value Cut-off control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43*

$$\text{Cut-off} = 0.43$$

### 9.3. Interpretation of Results

Samples are considered **POSITIVE** if the absorbance value is higher than 10% over the cut-off.

Samples with an absorbance value of 10% above or below the cut-off should not be considered as clearly positive or negative

→ **grey zone**

It is recommended to repeat the test again 2 - 4 weeks later with a fresh sample. If results in the second test are again in the grey zone the sample has to be considered **NEGATIVE**.

Samples are considered **NEGATIVE** if the absorbance value is lower than 10% below the cut-off.

#### 9.3.1. Results in NovaTec Units

$$\frac{\text{Patient (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec-Units} = \text{NTU}]$$

*Example:  $\frac{2.256 \times 10}{0.43} = 52 \text{ NTU (NovaTec Units)}$*

Cut-off:	10	NTU
Grey zone:	9-11	NTU
Negative:	<9	NTU
Positive:	>11	NTU

## 10. SPECIFIC PERFORMANCE CHARACTERISTICS

---

### 10.1. Precision

<b>Interassay</b>	<b>n</b>	<b>Mean</b>	<b>Cv (%)</b>
Pos. Serum	24	1.07	5.1
<b>Intraassay</b>	<b>n</b>	<b>Mean</b>	<b>Cv (%)</b>
Pos. Serum	8	1.12	2.6

### 10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is > 90 %.

### 10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 83,3 %.

### 10.4. Interferences

Interferences with hemolytic, lipemic or icteric sera are not observed up to a concentration of 10 mg/ml hemoglobin, 5 mg/ml triglycerides and 0.2 mg/ml bilirubin.

**Note:** The results refer to the groups of samples investigated; these are not guaranteed specifications

## 11. LIMITATIONS OF THE PROCEDURE

---

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values. Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

In immunocompromised patients and newborns serological data only have restricted value.

## 12. PRECAUTIONS AND WARNINGS

---

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The NovaLisa™ ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

WARNING:	In the used concentration Bronidox L has hardly any toxicological risk upon contact with skin and mucous membranes!
----------	---

WARNING:	Sulphuric acid irritates eyes and skin. Keep out of the reach of children. Upon contact with the eyes, rinse thoroughly with water and consult a doctor!
----------	--

### 12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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
## **BIBLIOGRAPHY / LITERATUR / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFÍA**

---

American Academy of Pediatrics. Parainfluenza Viral Infections. In: Peter G, ed. *1997 Red Book: Report of the Committee on Infectious Diseases*. 24th ed. Elk Grove Village, IL: American Academy of Pediatrics; 1997: 379.

Collins PL, Chanock RM, McIntosh K. Parainfluenza viruses. In: Fields BN, Knipe DM, Howley PM, eds. *Fields Virology*. 3rd ed. Philadelphia: Lippincott-Raven; 1995: 1205-41.

Glezen WP, Denny FW. Parainfluenza Viruses In: Evans A, Kaslow R, eds. *Viral Infections in Humans: epidemiology and control*. 4th ed. New York: Plenum; 1997:551-67.

Symbols Key/ Symbolschlüssel/ Explication des symboles / Legenda / Símbolos	
	Manufactured by / Hergestellt von/ Fabriqué par/ Prodotto da/ Fabricado por
<b>IVD</b>	In Vitro Diagnostic Medical Device/ In Vitro Diagnosticum/ Dispositif médical de diagnostic <i>in vitro</i> / Diganostico <i>in vitro</i> / Producto para diagnóstico In vitro
<b>LOT</b>	Lot Number/ Chargenbezeichnung/ Numéro de lot/ Lotto/ Número de lote
	Expiration Date/ Verfallsdatum/ Date de péremption/ Scadenza/ Fecha de caducidad
	Storage Temperature/ Lagertemperatur/ Température de conservation/ Temperatura di conservazione / Temperatura de almacenamiento
<b>CE</b>	CE Mark/ CE-Zeichen/ Marquage CE / Marchio CE/ Marca CE
<b>[REF]</b>	Catalogue Number/ Katalog Nummer/ Référence du catalogue/ Numero di codice/ Número de Catálogo
	Consult Instructions for Use/ Gebrauchsanweisung beachten/ Consulter la notice d'utilisation/ Consultare le istruzioni/ Consulte las Instrucciones de Uso
<b>MTP</b>	Microplate/ Mikrotiterplatte/ Microplaque/ Micropiastra/ Microplaca
<b>CONJ</b>	Conjugate/ Konjugat/ Conjugué/ Coniugato/ Conjugado
<b>CONTROL -</b>	Control serum, negative/ Kontrollserum, negative/ Sérum de contrôle négatif/ siero di controllo, negativo /Suero control negativo
<b>CONTROL +</b>	Control serum, positive/ Kontrollserum, positiv/ Sérum de contrôle positif/ siero di controllo, positivo/ Suero de control positivo
<b>CUT OFF</b>	Cut off control serum/ Cut off Kontrollserum/ Sérum de contrôle du cut-off/ siero di controllo, cut-off/ Suero control Cut-off
<b>DIL A</b>	Sample diluent buffer IgA/ IgA-Probenverdünnungspuffer/ Tampon diluant pour échantillon IgA/ soluzione tampone per i campioni IgA/ solución tampón para muestras IgA
<b>SOLN STOP</b>	Stop solution/ Stopplösung/ Solution d'arrêt/Soluzione bloccante
<b>SUB TMB</b>	TMB Substrate solution/ TMB-Substratlösung/ Substrat TMB/ soluzione substrato TMB/ solción substrato TMB
<b>WASHBUF 20x</b>	Washing solution 20x concentrated/ Waschlösung 20x konzentriert/ Solution de lavage concentré 20 x/ soluzione di lavaggio concentrazione x20/ solución de lavado concentrado x20
	Contains sufficient for "n" tests/ Ausreichend für "n" Tests/ Contenu suffisant pour "n" tests/ Contenido suficiente per "n" saggi/ Contenido suficiente para "n" tests

# SCHEME OF THE ASSAY

Parainfluenza Virus 1,2,3 IgA-ELISA

## Test Preparation

Prepare reagents and samples as described.  
Establish the distribution and identification plan for all specimens and controls on the result sheet supplied in the kit.  
Select the required number of microtiter strips or wells and insert them into the holder.

## Assay Procedure

	Substrate blank (e.g. A1)	Negative control	Positive control	Cut-off control	Sample (diluted 1+100)
Negative control	-	100µl	-	-	-
Positive control	-	-	100µl	-	-
Cut-off control	-	-	-	100µl	-
Sample (diluted 1+100)	-	-	-	-	100µl
Cover wells with foil supplied in the kit <b>Incubate for 1 h at 37°C</b> Wash each well three times with 300µl of washing solution					
Conjugate	-	100µl	100µl	100µl	100µl
Cover wells with foil supplied in the kit <b>Incubate for 30 min at room temperature</b> Wash each well three times with 300µl of washing solution					
TMB Substrate	100µl	100µl	100µl	100µl	100µl
<b>Incubate for exactly 15 min at room temperature in the dark</b>					
Stop Solution	100µl	100µl	100µl	100µl	100µl
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

## NovaTec Immundiagnostica GmbH

### Technologie & Waldpark

Waldstr. 23 A6  
D-63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760 Fax: +49 (0) 6074-487629

Email : [info@NovaTec-ID.com](mailto:info@NovaTec-ID.com)

Internet: [www.NovaTec-ID.com](http://www.NovaTec-ID.com)

PAIA0360engl,dt,fr,it,es01092011-CS